User Guide — SampleSheet Generator (WES, Pathology)

This app creates the Illumina `SampleSheet.csv` for a NovaSeq X Plus WES run. You enter each sample and choose its index; the app builds and checks the sheet for you and lets you export it.

If the app is not yet installed on your computer, see `docs/INSTALL.md` (installation is a one-time setup).

1. Start the app

1. Open the project folder.
2. Double-click `samplesheet-generator-pathology.Rproj`. RStudio opens.
3. In RStudio, open `app.R` (top-left file list) and click **Run App** (top-right of the editor).
4. The app opens in a window.

RStudio may briefly show a message about `renv` when the project opens. This is normal. As long as the app starts when you click **Run App**, everything is working.

2. Create a sample sheet

Enter the run name (optional)

- Type the **RunName** in the top field if you know it.
- If you don't know it yet, leave it empty — the sheet will contain `NA` in place of the run name, and the file will be named `NA-samplesheet.csv`. You (or Kai) can fill the run name in later.

Add a sample

For each sample:

1. Type the **Sample_ID** (for example `1234-26_3-N`).
2. In the **Index** field, start typing to search — you can search by the index name (e.g. `74`) or by sequence. Pick the correct index.
3. Click **Probe hinzufügen** (Add sample).

The sample appears in the table on the right, with its i7/i5 sequences filled in automatically. An index you have already used disappears from the list, so you cannot accidentally assign it twice.

The "**IndexUDP (nur Vorschau)**" column is shown only to help you check your work — it is **not** written into the sample sheet.

Remove a sample

Click the row in the table, then click **Ausgewählte Zeile entfernen** (Remove selected row). The index becomes available again.

3. Read the check ("Prüfung")

The app checks your samples continuously. You will see one of three states:

- **Green** — **"Keine Probleme gefunden."** Everything is fine; you can export.
 - **Yellow** — **Warnungen.** Something is unusual (for example a Sample_ID that does not match the expected pattern), but the sheet is still valid. A green note confirms: **"Der Lauf kann trotz Warnungen exportiert werden."** You can export.
 - **Red** — **Fehler.** There is a real problem (for example a duplicate Sample_ID, or an invalid character). **Export is blocked** until you fix it. Each message tells you which row and what is wrong.
-

4. Export the sample sheet

1. When the check shows green or yellow (not red), the **SampleSheet exportieren** button is active.
2. Click it. Choose where to save.
3. The file is saved as `<RunName>-samplesheet.csv` (or `NA-samplesheet.csv` if you left the run name empty).

If the button is grey, export is not possible yet — either you have no samples, or there are errors to fix (see the check panel).

5. Important note on accuracy

The app checks the **format** of each Sample_ID, but it **cannot** know whether you typed the correct number. A Sample_ID that is valid but wrong (for example `_1` typed instead of `_3`) will not be caught. **Always double-check the Sample_IDs against your source before exporting.** The review table is there to help you do this.

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